

# **REMOTE EVALUATION CHECK LIST**

# **IMPORTANT: PLEASE FOLLOW THESE INSTRUCTIONS TO ENSURE WE SUPPLY THE CORRECT BIDET.**


Section 1: Client Infor	mation											
Client Na	ame	Title/Pror			oun Phon		Number					
Postal Add	lress				'							
Suk	ourb	Postal Code				State Clier			nt Height/Weight			
Е	mail											
Section 2: Funding Inf	ormation											
Fund Gr	oup DVA	NDIS	0	Home Care P	ackage 🔘	Othe	er 🔘					
DVA File nun	nber		N	IDIS Numbei	-				DOB			
Plan Start I	Oate		Pla	n Finish Date	•			Funding	g Group			
E	mail											
Section 3: Prescribing	Therapist											
Occupational Thera	Phone Number			Number								
Е	mail											
Suk	ourb			Po	ostal Code					State		
Section 4: Toilet Meas	urements				NB: Our	bidet	s will n	ot fit on s	quare to	oilets		
PLEASE TA	KE PHOT	OS OF TOILI	ET (	O <sup>-</sup>		A					мм	
	<u> </u>	S	ee pag	e 2		<b>A1</b>					MM	
-		_				В					MM	
		—— A — — А 1 —		<b>→</b>		С					MM	
	ì	A.I				D					MM	
CISTERN	В					Е					MM	
	/ P-J	(← C −		)_ <sub>1</sub>	Where	is the v	water ta	p located?	OI	n Wall	In (	Cistern
TAP					Acces	s toile	t seat be	olts from	, Ur	nderneath	No	Access
E					Is there	access	to a po	wer point?	Y	es	No	0
	ls	there t	tank wa	ter in use?	Y	es*	No	O				
					*A photo o	of the wa	ater <b>tank</b>	<b>c pump</b> is re	quired.			
Please check with Cou In Australia, the plumbir					RPZ N	<u>//UST</u>	BEI	NSTALL.	E <mark>D.</mark> Ce	rtain ar	eas ac	cept
plumber and we do advi	80mm raiser, please check with local council.											
Must be installed	in compliance	with AS 3500.1						00 Series i	_			e
								nstallatio <mark>ther ba</mark>				
Section 5: ORDE	R FORM. Ple	ase tick if req	uired.		NONES							
Raisers	50mm	80mm			Arms	Ba	riatric	3in	1 🔲	Ende	eavour	
Pee Guard	Over	toilet chair		BA08 Tam	per Proof	Holde	er	В	A08/13	3 way c	ontrol	

For Aerolet Toilet Lift, please refer to the **Aerolet Evaluation form.** 




## PHOTO GUIDE - Below are examples of photos needed









#### Fitting measurements needed

- (A) Front (centre) of cistern to front (centre) of bowl
- (B) Centre to centre of bolt holes
- (C) Back (inside centre) to front (inside centre) of bowl
- (D) Bolt holes to front of bowl
- (E) Measure from the floor to top of toilet bowl

(A.1) If toilet bowl is not flat all the way to the cistern (lip) measure from the lip to front (centre) of bowl.

### **INSTALLATION & PLUMBING - AS 3500.1**

IF no 80m m raiser is being installed THEN an 'RPZ Valve' MUST be installed.

## Plumbing Regulations as per AUS/NZ standards

A toilet seat douche or bidet must comply with AS/NZS 3500 Part 1 Water Services Clause 12.2.3.

A douche outlet installed at least 25mm above the overflow level of the pan, at any position, does not require a backflow prevention device. If the douche outlet is less than 25mm above the overflow level of the pan, at any position, a 'high hazard' backflow prevention device must be installed. The device must comply with AS/NZS 2845.1 or AS 2845.2.

If InteliClean Air 9000 Series is prescribed, it meets the standard required installation AS2845.1 AS2845.2 & ASNZ3500. No further backflow is required.

An 80mm spacer may be sufficient. Please check with your local council/water regulator as this varies between council to council and state to state.

### SUMBIT FORM to OT@THEBIDETSHOP.COM.AU

If the client is DVA please submit this form and photos with the DOF and ATI.

Submit via email to ot@thebidetshop.com.au (or fax 07 5591 7744).

For assistance, please phone 1300 243 387 extension 3 (OT).